

Tech Note: Understanding a Certificate of Calibration

Introduction

This tech note describes information in the various sections of a Coherent certificate of calibration.


27650 SW 95th A Wilsonville, OR 9 Phone (800) 343-	ve. 7070			te of Ca 7025:2017 A		3/1/	CHILIT A C	ANAB Rational Accreditation Board C R E D I T E D SOINC 17025 BRATION LABORATORY
Date: 29	Sep 2023					Certifica	ation Number: 2309	29135322
Part Number: 001	12-1180					Tem	perature (°C): 22.7	
Description: PM	200F-50					Relative	Humidity (%): 35.7	
Serial Number: 000	08A05						Procedure: QI-19	9.70 RevGG
Instrument Co	ondition As F	Receive	d					
Wavelength	Responsivity	Unce	onsivity rtainty =2)	Laser Power	Measured Powe	Measured Power Uncertainty* (k=2)	Laser Power Tolerance Limit	Status
514 nm	6.206E-4 V/W	±1.0 %		1.004W	1.007W	±1.4 %	990.1mW - 1.018W	In Tolerance
Instrument Co Wavelength	Responsivity 6.206E-4 VW	Respo	onsivity ertainty =2)	Laser Power	Measured Power	Measured Power Uncertainty* (k-2)	Laser Power Tolerance Limit 990.1mW -	Status In Tolerance
	Power Uncertainty Standards	" includes	the uncer	tainty of the meter u	sed for this meas		ation Due	
MOLECTRON P	M10		0011W00	0	I	Dec 2023		
Coherent LabM	ax TOP		1781T12	!	/	Apr 2024		
HP 34401A			0124T97	,	1	Nov 2023		
Comments:	C BLACK			Т	est Technician			
The results listed o written approval by using standards tra	rval begins when Date". Contact Contac	ustomer S certificate a tion certific rnational S	ervice for i applies to cates with system of t	recommended callt only the item listed out signatures are r Units (SI) via a Natio	pration intervals for above and shall r not valid. Coheren onal Metrology Ins	or Coherent products not be reproduced of it hereby certifies that stitute (NIST, PTB, et	s. ther than in full, with at the above item ha c.), that are signator	s been calibrated les to the CIPM
a coverage factor of The acceptance cri	f k=2. teria for the specif	ication(s) i	n this cert		er ILAC-G8:09/20	119, sec. 4.2.1 - Bina		onfidence level using nple Acceptance



Coherent ANAB Accreditation and Calibration Date

The header section displays the Coherent address and standards to which Coherent is accredited, including the ISO/IEC 17025:2017 standard. The ANAB logo verifies that the Coherent product, listed on the calibration certificate, is accredited and certified. If the ANAB logo is missing signifies a Coherent product outside the Coherent ANAB scope of accreditation.

The date, stated on the top of the certificate, is the date of calibration.



Instrument Condition As Received

This section identifies the product's As Received condition as returned to Coherent for calibration. The As Received lists the operating conditions and status under which the unit is tested upon receipt.


A blank As Received identifies the calibration certificate as the initial production calibration certificate or that the unit as received state was unable to provide a valid as received condition. This happens when the element is damaged or broken beyond a reasonable measurement collection.

Wavelength	Responsivity	Responsivity Uncertainty (k=2)	Laser Power	Measured Power	Measured Power Uncertainty* (k=2)	Laser Power Tolerance Limit	Status



Instrument Condition As Shipped

The As Shipped section describes the final calibration condition of the unit when shipped to the customer.




- 1. The Laser Power/Energy is for this specific measurement, as measured by the Working Standard.
- 2. The Measured Power/Energy is the measurement made by the Unit Under Test.

measurement, by the Unit Under Test, is within these limits it is in tolerance.

- 3. Uncertainty values include both the meter and sensor, added up using the sum of squares method. Therefore, the uncertainty is slightly higher than the uncertainty for the sensor Responsivity (Rv) value only. (See the Note on the Certificate that the "Measured Power Uncertainty" includes the uncertainty of the meter used for this measurement.)
- 4. The Tolerance Limit is based upon applying the uncertainty of the measurement to the laser Power/Energy, as measured by the Working Standard.
 Upper Tolerance Limit = Laser Power + (Laser Power * Uncertainty)
 Upper Tolerance Limit = Laser Power (Laser Power * Uncertainty)
 This lists the allowable range of Measured Power/Energy or the "tolerance limits". If the
- 5. The Status column compares the Measured Power/Energy to the Tolerance Limits and reports if the Measured Power/Energy is In Tolerance or Out of Tolerance.

Equipment Used for Calibration

This section lists the Coherent equipment used to perform the calibration and against which the calibration is measured for standard traceability. The standard equipment Calibration Due date is listed and understood to be on the last day of the month.





Calibration Decision Rule

The Calibration Interval Start Date and Due Date are left blank for the customer to identify the recalibration requirements after the product is placed in service.

The decision rule explains the calibration acceptance criteria for the Tolerance specifications.

Coherent reports the Uncertainty at the k=2 confidence level. This is also called an Expanded Uncertainty, defined as two Standard deviations. This equates to 95% confidence interval.

alibration Interval Start Date: Due

Due Date:

The calibration interval begins when the equipment is placed into service. The "Due Date" may be established (by the customer) by adding the calibration interval to the "Start Date". Contact Customer Service for recommended calibration intervals for Coherent products.

The results listed on this calibration certificate applies to only the item listed above and shall not be reproduced other than in full, without the specific written approval by Coherent. Calibration certificates without signatures are not valid. Coherent hereby certifies that the above item has been calibrated using standards traceable to the International System of Units (SI) via a National Metrology Institute (NIST, PTB, etc.), that are signatories to the CIPM Mutual Recognition Agreement. The reported uncertainties are expressed as expanded uncertainty values at an approximately 95% confidence level using a coverage factor of fe2.

The acceptance criteria for the specification(s) in this certificate are defined per ILAC-G8.09/2019, sec. 4.2.1 - Binary Statement for Simple Acceptance Rule: In Tolerance - Compliant with specification(s); Out of Tolerance - Not compliant with specification(s).

Page 1 of 1

Contact Coherent

For assistance or additional information, contact Coherent Technical Support as follows:

- Contact your local Coherent Service Representative (or visit <u>www.Coherent.com</u> to view a list of contacts worldwide)
- Send an e-mail to: <u>LSMservice@Coherent.com</u>
- Call the Coherent Technical Support Hotline

Within the USA: 1-(800)-343-4912
 Outside of the USA: 1-(408)-764-4042

Additional sensor products information can be found on the Coherent website under Lasers/Laser Measurement Instruments or this address:

https://www.coherent.com/laser-power-energy-measurement

Sensor product current software downloads can be found on the Coherent website under Support/Downloads or this address:

https://www.coherent.com/resources?resourceType=Manual%2BSoftware%2BLMC%20Video%20Library

To arrange for warranty service or annual recalibration, contact your regional Coherent service center to obtain a Return Material Authorization (RMA) number. Use the shipping box and packaging materials you retained to safely transport the sensor back to the factory, and ship to this address:

Coherent, Inc. Attn: RMA # 27650 SW 95th Ave. Wilsonville, OR 97070

