


CABLE LENGTH 1220mm ±25mm [ 48" ±1"]





**CONNECTOR**

6-PIN MALE CONNECTOR  
MOLEX PART NUMBER 43025-0600

6X CONTACT PIN,  
MOLEX PART NUMBER 43030-009

**SCHEMATIC**

PINS, 1,2,3, GROUND  
PINS 4,5,6 AT 12 VDC



189.0mm  
[ 7.44in ]

89.4mm  
[ 3.52in ]

47.1mm  
[ 1.85in ]

ON/OFF SWITCH

IEC 320-C14 INPUT RECEPTACLE  
UNIVERSAL 90/264VAC CAPABLE INPUT


Title:  
**OBIS Six Laser Power Supply**



**COHERENT**  
5100 Patrick Henry Drive  
Santa Clara, CA 95054  
www.coherent.com

PROPRIETARY NOTICE:  
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF  
COHERENT, INC. IT IS NOT TO BE USED, REPRODUCED, OR DISCLOSED IN WHOLE  
OR IN PART TO ANYONE WITHOUT THE WRITTEN PERMISSION OF COHERENT, INC.

DRAWING NOT TO SCALE    Date: 01-09-12    Rev. 01



THIRD ANGLE  
PROJECTION