



CERTIFICATE



This is to certify that

CoAdna (Suzhou) Limited

No.4 Genway Factory, No 158, Qiming Road, Comprehensive Bonded Zone, Suzhou Industrial Park
215121, Suzhou, Jiangsu
China

has implemented and maintains a **Quality Management System**.

Scope:

Manufacturing of Wavelength Selective Switch for Optical Networks

Product Categories:

Product Category: 7.7.4 - Electromechanical
Assemblies

Specialty Option: Hardware

Exemptions: NA

Through an audit documented in a report, in accordance with QuEST Forum Code of Practice for TL 9000 Registrars and ISO/IEC 17021 Conformity Assessment - Requirements for bodies providing audit and certification of management systems, it was verified that the management system fulfills the requirements of the following standards:

TL 9000-H R6.3/R5.7

ISO 9001 : 2015

TL Id number	TL 6174
Certificate registration no.	50052711 TLR6
Date of certification	2021-11-11
Valid until	2024-11-10



DQS Inc.

Brad McGuire
Managing Director

Accredited Body: DQS Inc., 1500 McConnor Parkway, Suite 400, Schaumburg, IL 60173 USA
Administrative Office: DQS AP Ltd., 906-907, Waterfront Place Block E, No.31,Lane 168, Daduhe Road,
Putuo District, Shanghai, China, Post Code: 200062